


# Metagenomic Sample Analysis

## Standard Operating Procedures for fossil fuel environments




HYDROCARBON METAGENOMICS

## Contents

1	Sample Collection .....	4
1.1	Metadata.....	4
1.2	Sample Fractionation .....	4
2	DNA Extraction.....	5
2.1	DNA Sample Fractionation.....	7
3	Community Analysis.....	8
3.1	PCR amplification .....	8
	Procedure – PCR amplification and barcoding .....	9
	Materials and Reagents .....	9
3.2	16S DNA Sequencing.....	11
3.3	Data Analysis via Phoenix II.....	11
4	Total DNA Analysis .....	12
4.1	Total DNA Sequencing .....	12
4.2	Binning and Assembly .....	12
4.3	Annotation .....	12
4.4	Data Analysis .....	12
5	Data Storage.....	13
6	Appendix 1 .....	14

## Typical Metagenomics Work Flow



# 1 SAMPLE COLLECTION

---

The goal of effective sample collection is to minimize contamination while maintaining the character of the indigenous microbial community and obtaining sufficient amounts of biomass. The sample and resulting sequence data should be checked for contamination.

The sample procedure will be dependent on the nature of the natural environment.

- I. Water (Produced water, Tailings pond water, Pipeline water)
  - Sample should be collected in sterile Nalgene bottles filled to the top to avoid air exposure. When samples are processed immediately they can be stored in an anaerobic chamber containing 90% N<sub>2</sub> and 10% CO<sub>2</sub>. Samples to be kept for later metagenomic analysis should be stored at -80 °C. Tailings ponds often yield around 100 ng per ml of sample, while produced water yields are lower, at around 5 ng per ml.
  - Excess water from tailings samples are removed with centrifugation at 14,000 rpm
  - Produced water samples are filtered with a 0.2µm filter
  - Oil can be removed from samples by a short centrifugation at 14,000 rpm
- II. Drilling Cores (Coal, Oil Sands, etc)
  - Drilling cores should be rapidly frozen at the well site immediately following collection. Cores should remain frozen when in transportation. Samples should be stored at -80°C. Core can be subdivided using a sterile rock saw. Samples should be prepped for DNA extraction by aseptically removing the outside surfaces of the cores (to eliminate drilling contamination) leaving the interior of the core material for DNA analysis. The interior frozen core material is then ground with a mortar and pestle prior to DNA extraction. The yield of core material is often less than 1 ng of DNA per g of core material. Although freezing is ideal for metagenomics (DNA) analysis it is not appropriate for all downstream analysis, such as culturing. If other tests will be conducted on the core material, some of the core material should be left unfrozen.
- III. Biofilms (pipelines)
  - Pipeline sections should be stored in the anaerobic hood (90% N<sub>2</sub> and 10% CO<sub>2</sub>) until the pipeline biofilm can be scraped off. The scrapings should be processed immediately or stored at -80°C. Although the yield can vary widely, based on the composition of the scrapings, it can be as high as 500 ng per g.

## 1.1 METADATA

Effective sample collection also requires extensive metadata collection. Items that should be recorded include: Environmental conditions (temperature, depth, and location), Chemistry (pH, salinity, etc), and sampling methodology to name a few. More information on metadata collection can be found at the Genomics Standards Consortium's "Minimum Information about a (Meta)Genome Sequence" (MIMS) website: [www.gensc.org](http://www.gensc.org).

## 1.2 SAMPLE FRACTIONATION

As long as sufficient material is obtained, samples should be subdivided and some material stored at -80 °C for later use or for different analyses.

## 2 DNA EXTRACTION

---

The optimal DNA extraction procedure will depend on the sample being processed. For the majority of samples a bead-mill homogenization with SDS-Chloroform Treatment can be used for DNA extraction. The protocol is listed below, other protocols commonly used are listed in appendix 1.

### **DNA extraction from tailings ponds using FastDNA spin kit for soil**

#### **Overview**

The FastDNA® SPIN Kit for Soil quickly and efficiently isolates PCR-ready genomic DNA directly from soil type samples in less than 30 minutes. Designed for use with the FastPrep® Instruments from MP Biomedicals, microbial cells within soil type samples are easily lysed within 40 seconds. These benchtop devices use a unique, optimized motion to homogenize samples by multidirectional, simultaneous impaction with lysing matrix particles.

**Part I** describes the cell lysis step and the isolation of genomic DNA from soil type samples. Samples are placed into 2.0 ml tubes containing Lysing Matrix E, a mixture of ceramic and silica particles designed to efficiently lyse all soil organisms including historically difficult sources such as eubacterial spores and endospores, and gram positive bacteria. Homogenization in the FastPrep® Instrument with Lysing Matrix E takes place in the presence of MT Buffer and Sodium Phosphate Buffer, reagents carefully developed to protect and solubilize nucleic acids and proteins upon cell lysis. These reagents work together to allow extraction of genomic DNA with minimal RNA contamination. Following lysis, samples are centrifuged to pellet soil, cell debris and lysing matrix. DNA is purified from the supernatant with a silica-based GENCLEAN® procedure using SPIN filters. **Part II** describes the quality control step, where the resulting DNA is assessed for its quantity, quality, and purity. This would determine whether the DNA is appropriate for 454 pyrosequencing.

#### **Kit Components and User Supplied Materials**

##### **FastDNA® SPIN Kit for Soil Components**

Lysing Matrix E	50x 2.0 ml tubes
Sodium Phosphate Buffer	60 ml
MT Buffer	8 ml
PPS Solution	25 ml
Binding Matrix	66 ml
SPIN Modules	50 each
Catch Tubes	50 each
Concentrated SEWS-M	2 ml
DES	20ml
BBS Gel Loading Dye	200 µl
User manual	1 each
MSDS	1 each
Certificate of Analysis	1 each

### User Supplied Materials

FastPrep® Instrument  
Microcentrifuge that can freely spin 2.0 ml tubes  
Microcentrifuge tubes (2.0 ml and 1.5 ml)  
Clean 15 ml tubes for DNA binding  
Rotator or low-speed vortex

### Reference:

<http://images.www.mpbio.com/docs/fastprep/FastDNASpinKitforsoil.pdf>

### Part I: DNA Extraction

1. Add up to 500 mg of sample to a Lysing Matrix E tube.

If there is too much water in the sample, it is possible to remove the supernatant with centrifugation, which may lead increased DNA yields.

2. Add 978 µl Sodium Phosphate Buffer to sample in Lysing Matrix E tube.
3. Add 122 µl MT Buffer.
4. Homogenize in the FastPrep® Instrument for 40 seconds at a speed setting of 6.0.
5. Centrifuge at 14,000 x g for 5-10 minutes to pellet debris.  
NOTE: Extending centrifugation to 15 minutes can enhance elimination of excessive debris from large samples, or from cells with complex cell walls.
6. Transfer supernatant to a clean 2.0 ml microcentrifuge tube. Add 250 µl PPS (Protein Precipitation Solution) and mix by shaking the tube by hand 10 times.
7. Centrifuge at 14,000 x g for 5 minutes to pellet precipitate. Transfer supernatant to a clean 15 ml tube.  
NOTE: While a 2.0 ml microcentrifuge tube may be used at this step, better mixing and DNA binding will occur in a larger tube.
8. Resuspend Binding Matrix suspension and add 1.0 ml to supernatant in 15 ml tube.
9. Place on rotator or invert by hand for 2 minutes to allow binding of DNA. Place tube in a rack for 3 minutes to allow settling of silica matrix.
10. Remove and discard 500 µl of supernatant being careful to avoid settled Binding Matrix.
11. Resuspend Binding Matrix in the remaining amount of supernatant. Transfer approximately 600 µl of the mixture to a SPIN™ Filter and centrifuge at 14,000 x g for 1 minute. Empty the catch tube and add the remaining mixture to the SPIN™ Filter and centrifuge as before. Empty the catch tube again.
12. Add 500 µl prepared SEWS-M and gently resuspend the pellet using the force of the liquid from the pipet tip.  
NOTE: Ensure that ethanol has been added to the Concentrated SEWS-M.
13. Centrifuge at 14,000 x g for 1 minute. Empty the catch tube and replace.
14. Without any addition of liquid, centrifuge a second time at 14,000 x g for 2 minutes to “dry” the matrix of residual wash solution. Discard the catch tube and replace with a new, clean catch tube.
15. Air dry the SPIN™ Filter for 5 minutes at room temperature.
16. Gently resuspend Binding Matrix (above the SPIN filter) in 50-100 µl of DES (DNase/Pyrogen-Free Water).

NOTE: To avoid over-dilution of the purified DNA, use the smallest amount of DES required to resuspend Binding Matrix pellet.

NOTE: Yields may be increased by incubation for 5 minutes at 55°C in a heat block or water bath.

17. Centrifuge at 14,000 x g for 1 minute to bring eluted DNA into the clean catch tube. Discard the SPIN filter. DNA is now ready for PCR and other downstream applications. Store at -20°C for extended periods or 4°C until use.

### Part II: Quality Control

1. Quantify the DNA concentration via Quant-iT™ fluorometer.
2. Ensure DNA can be amplified with PCR primers.

## **2.1 DNA SAMPLE FRACTIONATION**

DNA should be subdivided so some DNA can be stored frozen at -80 °C for later use.


### 3 COMMUNITY ANALYSIS

#### 3.1 PCR AMPLIFICATION

Either primers that universally amplify the 16S gene can be used or primers targeting specific organisms or functional genes can be used. It is important to know the limitations (biases) of the primers you are using.

Since the PCR primers with barcodes and sequencing (FLX) adapters are quite long (40 nt), PCR amplification is performed in two stages. First the DNA is amplified with only the universal primers (no adapters) and the sequencing adapters are added during a second round of PCR.

**Workflow for two-step PCR amplification using 454-sequencing primers**



**Universal PCR PRIMERS** (See Roche technical bulletin [TCB No. 005-2009](#) for list of bar codes (MIDs))

- Original PCR primer no adaptor

1392R: acg ggc ggt gtg tRc

926wF: aaa ctY aaa Kga att gRc gg

- Original PCR primer with adaptor

_____ FLX Adapter _____	_____ bar code _____	_____ Primer _____
454T_RL 5'-CCATCTCATCCCTGCGTGCTCTCCGAC	TCAG ACGAGTGCGT	acg ggc ggt gtg tRc
454T-FB 5'-CCTATCCCCTGTGTGCCTTGCCAGTC		aaa ctY aaa Kga att gRc gg

- Updated PCR primer for Epsilon no adaptor

1392Re: acg ggc ggt gWg tRc

926wFe: aaa ctY aaa Kga atW gRc gg

- Updated PCR primer for Epsilon with adaptor

_____ FLX Adapter _____	_____ bar code _____	_____ Primer _____
454T_RL 5'-CCATCTCATCCCTGCGTGCTCTCCGAC	TCAG ACGAGTGCGT	acg ggc ggt gWg tRc
454T-FB 5'-CCTATCCCCTGTGTGCCTTGCCAGTC		aaa ctY aaa Kga atW gRc gg



## Procedure – PCR amplification and barcoding

- Establish unidirectional flow of reagents.
- Clean work surfaces and tools with 70% ethanol before and after each use

### Materials and Reagents

Starting materials

-Extracted DNA diluted to 10ng/ul per sample

<u>Materials/Reagents/Equipment</u>	<u>Vendor</u>	<u>Stock Number</u>
Taq polymerase 10x PCR buffer included	Amersham Pharmacia Biotech	27-0799-62
BSA	New England Biolab	B9001S
Nuclease Free Water (500mL)	Ambion	9930
10mM dNTP Mix (200ul, 1ml)	MBI Fermentas	R0191, R0192
100bp DNA ladder Plus (0.5ug/ul, 50ug)	MBI Fermentas	SM0323
Minelute Gel Extraction Kit (50,250)	Qiagen	28604, 28606
TE		
Seakem GTG agarose		
Ethanol		
<u>Equipment</u>		
MJ Thermocycler (96-well alpha)	MJ Research	ALS 1296
Microcentrifuge (for 1.5 tubes)	ISC Bioscience	C-1201
Microcentrifuge (for 8-strip tubes)	ISC Bioscience	C-1202
Vortex Genie	ISC Bioscience	S-7350-1
Pipettes (P2, P10, P20, P200, P1000)	Rainin	RL-2, 10, 20, 200, 1000
0.2mL tube rack	Fisher	05-541-55
8-tip multichannel pipettor	Rainin	L8-200
<u>Disposables</u>		
8-strip tubes and caps (0.2mL)	ISC Bioscience	T-3014-1
0.2 mL tubes (with cap)	ISC Bioscience	C-3328-1
Tips (P2, P10, P20)	Rainin	RT-L10F
Tips (P200)	Rainin	RT-L200F
Tips (P1000)	Rainin	RT-L1000F
Tips (P1000)	Rainin	RT-L1000F

### 1. Prepare oligonucleotides

- 1.1. Order 454 primer sequences, 16S sequences, and barcodes.
  - a. Scale oligos to 25nmole, purified using standard desalting
- 1.2. Upon arrival, briefly spin down tubes containing primers in a minicentrifuge as the contents may have become loose during shipment
  - a. Resuspend each primer in TE to a stock concentration of 100uM

### 2. Prepare PCR cocktail.

- 2.1 For each sample to be amplified, prepare mix listed below. To prevent depletion of mix due to pipeting error, add 3 additional reactions to your total.

- a. The following is a list of reagents added **per reaction**. PCR master mix should be prepared as follows;
  - 13.9 uL water
  - 2.0 uL 10x PCR buffer
  - 0.4 uL 10mM dNTPs
  - 0.6 uL 10mg/ml BSA
  - 1.0ul 1uM Forward primer
  - 0.1 uL Taq
- 2.2 Vortex to mix, spin down
- 2.3 Dispense 18ul aliquots of the master mix into 0.2mL microtubes
- 2.4 Add the following in **triplicate** tubes/wells
  - 1ul of the 1uM Reverse primer
  - 1ul 10ng/ul DNA
- 2.5 Vortex to mix, spin down, cap.
- 3 Enter the following PCR program into Thermal cycler. See Thermal cycler manual for instructions on entering new PCR programs.
  - 95C, 3 min
  - 25 cycles of: 95C 30s, 55C 45s, 72C 90s \*
  - 72C 10min
  - 4C ∞
  - \*second round of amplification with adapters only uses 10 cycles*
- 4 Confirm successful amplification by running reactions out on a 1% agarose gel
- 4.1 Preparing Agarose Gel and run buffer
  - a. Make a 1.0% GTG SeaKem agarose TAE mini gel
  - b. Stir 200ml 1X TAE with 2.0g agarose and microwave until completely melted, stir on stirplate until slightly warm then add 10ul EtBr.
  - c. Pour 1X TAE buffer to fill line of gel tank.
- 4.2 Running Agarose Gel
  - a. Submerge solidified gel in TAE run buffer in tank, and remove comb.
  - b. Add 5-6ul of 5X loading dye to each sample, vortex to mix.
  - c. Load Samples in large wells and add a 100bp ladder
  - d. Run gel at 120V until band size is easily visible ~30-40 minutes
  - e. Image gel using dark reader settings on Fluor-S multi imager.
- 5 Purify amplification (PCR) product by QIAquick PCR Purification (QIAGEN)/ other similar kit
  - a. Add 1:250 volume of pH indicator I to Buffer PB, depending on the amount to be used (e.g. 16 µl pH indicator I to 4 ml Buffer PB). Yellow color of Buffer PBI indicates pH ≤ 7.5.
  - b. Next, add 5 volumes of Buffer PBI to 1 volume of the PCR sample and mix (e.g. add 500 µl of Buffer PBI to 100 µl PCR sample).
  - c. If the color of the mixture is (pale) orange or violet, add 10 µl of 3 M sodium acetate, pH 5.0 (need to be prepared, not included in the kit) and mix. The color of the mixture will turn to yellow.
  - d. Place a QIAquick spin column in a provided 2 ml collection tube.
  - e. To bind DNA, apply the sample to the QIAquick column and centrifuge for 30-60 s.
  - f. Discard flow-through. Place the QIAquick column back into the same tube.
  - g. To wash, add 750 µl Buffer PE to the QIAquick spin column and centrifuge for 30-60 s.

- h. Discard flow-through and place the QIAquick column back in the same tube. Centrifuge the column for an additional 1 min.
- i. Place QIAquick column in a clean 1.5 ml microcentrifuge tube.
- j. To elute DNA, add 30 µl - 50 µl Buffer EB (10 mM Tris-Cl, pH 8.5) or water (pH 7.0-8.5) to the center of the QIAquick membrane and centrifuge the column for 1 min.

*Addition of the elution buffer can vary depending on the amount of the PCR product obtained.*

6. Quantify each purified sample


[For additional details, please refer QIAGEN/QIAquick Spin Handbook (March 2008).]

### 3.2 16S DNA SEQUENCING

In order to obtain read lengths long enough to identify taxa to the genus level, Roche's 454 flex titanium platform is used. Additionally the length of the primers used in the PCR amplification are optimized for the 500 bp reads of the 454 flex.

### 3.3 DATA ANALYSIS VIA PHOENIX II

Once the 16S sequence data is available the data needs to go through quality control to eliminate the poor or chimeric sequences. The quality controlled data is then aligned and clustered so diversity statistics and taxonomic identifications can be made. The following figure demonstrates a typical analysis workflow.



There are several tools that can do this analysis. These include Phoenix 2 (<http://hmp.ucalgary.ca/phoenix/>) and QIIME (<http://qiime.org/>).

## 4 TOTAL DNA ANALYSIS

---

Total DNA analysis avoids the PCR amplification stage and is sent directly for sequencing.

### 4.1 TOTAL DNA SEQUENCING

Currently two sequencing platforms are used for total DNA sequencing

1) Roche 454 flx+

- the ~800 bp reads provide more reliable annotation for unassembled metagenomes.
- ½ run of 454 sequencing typically provides 800,000 reads (640,000,000 bps)

2) Illumina HiSeq 2000

- All though the reads from the HiSeq 2000 are much shorter ~150 to 200 bps, paired end runs are possible and the sequencing depth is much higher with 1.5 billion reads per 2 lanes (~40,000,000,000 bps)
- the greater coverage allows for the identification of rarer genes and species but it makes assembly and annotation more difficult.

### 4.2 BINNING AND ASSEMBLY

Binning and assembling large metagenomics datasets is complex and requires a dedicated and sophisticated computing systems. However metagenomic datasets can be and are often analyzed without assembling the dataset. Besides being computationally complex, assembling the data can create chimeric sequencing. Unassembled data (gene centric analysis) can be conducted using the IMG/MER and MG-Rast platforms described in 4.3 & 4.4. The project's binning and assembly SOP is below.

Project SOP coming soon

### 4.3 ANNOTATION

Project SOP coming soon.

Although Online tools like MG-Rast (<http://metagenomics.anl.gov/>) and IMG/MER (<https://img.jgi.doe.gov/mer/>) will not assemble metagenomics data for end users, they will annotate your datasets. Programs like MLTreeMap (<http://www.mltreemap.org/>) will use conserved genes provide taxonomic assignment to the microbes in the metagenome sample.

### 4.4 DATA ANALYSIS

Both MG-Rast and IMG/MER provide a wide range of tools for metagenome data analysis and comparison tools.

## 5 DATA STORAGE

---

Although raw data should be stored on site and remotely backed up, online databases like the NCBI's SRA will accept both 16S amplicon data and total DNA metagenomes datasets. Additionally, both MG-rast and IMG/MER will store the quality controlled and annotate reads from metagenome samples.

## 6 APPENDIX 1

---

Alternative DNA extraction protocols:

### **DNA extraction using chloroform, isopropanol and bead beating**

#### **Materials:**

- Solutions (all sterile filtered):
  - 50mM phosphate buffer pH 7.0
  - lysis buffer: 10% SDS, 100mM NaCl, 500mM Tris pH 8.0
  - 24:1 chloroform:isoamyl alcohol (make reasonably fresh)
  - 7M ammonium acetate
  - 100% isopropanol
  - MilliQ-filtered water
- Equipment:
  - 1.5ml eppendorff tubes
  - 2ml screw-cap tubes containing 0.5g (ea.) of 0.5 and 10 micron diameter zirconium beads
  - pipettes and tips
  - bead beater
  - microcentrifuge

#### **Method:**

- 1) In a 2ml screw-cap tube containing 0.5g (ea.) of 0.5 and 10 micron diameter zirconium beads add 300µl 50mM phosphate buffer pH 7.0
- 2) Add a healthy loop full of cells
- 3) Add 300µl lysis buffer and 300µl 24:1 chloroform:isoamyl alcohol
- 4) Shake on a bead beater for 40 seconds at 5000rpm
- 5) Centrifuge at 13,000rpm for 5 minutes
- 6) Remove the supernatant (~650µl) into a fresh 1.5ml eppendorff tube
- 7) Add 7M ammonium acetate to a final concentration of 2.5M (~360µl) and gently mix
- 8) Centrifuge at 13,000rpm for 5 minutes
- 9) Remove the supernatant (~580µl) into a fresh 1.5ml eppendorff tube
- 10) Add 1 volume of 100% isopropanol
- 11) Incubate at -20°C for at least 30 minutes
- 12) Centrifuge at 13,000rpm for 10 minutes
- 13) Remove supernatants and dry pellets
- 14) Resuspend in 40µl MilliQ water

#### **Notes:**

- Procedure works well done in batches of 12 reactions (capacity of the bead-beater)
- Total procedure up to incubation at -20°C takes approximately 45 minutes
- Extracts can be incubated at -20°C at least overnight (indefinitely?)

# **DNA EXTRACTION for HIGH MOLECULAR WEIGHT DNA USING 0.22 µM STERIVEX FILTERS**

Jody J. Wright, Elena Zaikova, David A. Walsh and Steven J. Hallam

Department of Microbiology & Immunology

University of British Columbia

**Corresponding author:** Steven J. Hallam

**Keywords:** microbiology, seawater, DNA extraction

**Short Abstract:** We describe a method for extraction of high molecular weight genomic DNA from planktonic biomass concentrated on 0.22 µm Sterivex filters.

**Long Abstract:** This method is used to extract high molecular weight genomic DNA from planktonic biomass concentrated on 0.22 µM Sterivex filters that have been treated with storage/lysis buffer and archived at -80°C. The protocol begins with two one-hour incubation steps to liberate DNA from cells and remove RNA. Part II involves a series of Phenol:Chloroform and Chloroform extractions followed by centrifugation to remove proteins and cell membrane components, collection of the aqueous DNA extract, and several buffer exchange steps to wash and concentrate the extract. Part III describes a method to check for the quality and presence of DNA in the extract and recommendations for both long and short-term storage of the resulting samples. The total time required for this protocol depends on the number of samples to be extracted. It is recommended to work with less than 15 samples at one time to avoid confusion and cut down protocol time. For 10-15 samples and assuming the proper centrifugation equipment is available, this protocol should take about 8 hours, with an additional overnight step to run an agarose gel to visualize and quantify the DNA extracts. Make sure you have the hybridization ovens set to temperature at the outset of the process.

## **Protocol**

### **Part 1: Incubation**

- 1) Thaw filters on ice.
- 2) Add 100 µl lysozyme (125 mg in 1000 µl TE) and 20µL RNase A (10 µg/ml) to each filter. Reseal the top with Parafilm. Incubate in a rotating incubator at 37°C for 1 hour.
- 3) Add 100µl Proteinase K and 100 µl 20% SDS to each filter. Reseal using Parafilm. Incubate at 55°C for 1 – 2 hours in a rotating incubator.
- 4) Remove lysate from sterivex filter using a 5cc syringe into a 15 ml falcon tube. Add 1 ml lysis buffer to rinse the fillter out. Pool this with the lysate in the 15 ml tube.

### **Part II: Extraction & Centrifugation**

- 5) Add an equal volume (about 3ml) of Phenol:Chloroform:IAA (25:24:1), pH 8.0 to the lysate tube. Vortex for 10 seconds. Spin at 2500 g for 5 minutes or until aqueous layer is clear (use the J.S.-5.3 rotor, using the 50 ml falcon tubes in which the filter was stored as an adapter). Transfer aqueous layer into a new 15 ml falcon tube.

- 6) Add an equal volume (approx 3mL) of Chloroform:IAA (24:1) to the tube containing the aqueous layer. Vortex for 10 seconds. Spin at 2500 g for 5 minutes or until aqueous layer is clear with no debris. Transfer aqueous layer into a new, labeled falcon tube. Add 1ml of TE buffer (pH 8.0) to this tube.
- 7) Transfer contents of the falcon tube for step 6) to an Amicon Ultra centrifuge tube. Spin at 3500 g for 10 minutes. Check to make sure there is less than 1ml liquid left in Amicon at the end of this (if not, fill up with TE and spin again). Remove filtrate to another falcon tube and save in fridge until DNA has been recovered on gel.
- 8) Add 2 ml TE buffer to Amicon and spin at 3500 g for 6 minutes. Remove filtrate.
- 9) Repeat 8) twice more (total of 3 washes with 2 ml TE). For the last wash, spin until 200 – 500 µl remain in Amicon. Note the final volume and transfer to a labeled 1.5 µl eppendorf tube.

### **Part III: Storage & Confirmation of DNA Presence**

- 10) Aliquot 70µl from final volume into a 1.5µl eppendorf tube to use as working stock and place in -20°C freezer. Place the rest of the DNA into -80°C freezer.
- 11) Set up an overnight agarose gel (0.8%) loading 1-2 lanes of 10 µl 1kb+ ladder and 3 lanes of 50ng/µl λHindIII ladder (2µl, 5µl, and 10µl) to use as size and intensity standards. Load 5 µl per lane of DNA extract sample. Run gel at 15 volts for approximately 16 hours to and check image in UV gel dock for presence of DNA in extracts.

**Representative Results:** – When this protocol is done correctly, you should see a gel image similar to Figure 1. Actual DNA concentration of extracts will vary depending on the source of the sample.

**Discussion:** Depending on how many samples are to be processed, this can be a time-intensive procedure due to the two one-hour incubation steps, and the repeated washes and centrifugation steps. It is best to plan a whole day for this procedure to leave plenty of time. If there are problems getting the final extract volume in the Amicon tube to reduce down to 200-500µl, try doing additional TE washes and centrifugations (repeat step 9) until the appropriate volume is leftover. Also, if the extract smells acutely of chloroform, it is best to do additional washes until the smell lessens to make sure all of the chloroform is removed.


### **Table of specific reagents and equipment:**

<b>Name of the reagent</b>	<b>Company</b>	<b>Catalogue number</b>	<b>Comments (optional)</b>
Lysozyme	Sigma	L6876-5G	
RNase A	Fisher	BP2539-100	
Proteinase K	Qiagen	19133	
Lysis Buffer			0.75 M sucrose, 40 mM EDTA, 50 mM Tris (pH 8.3)



Parafilm	Fisher	13-374-10	
20% SDS			
5cc syringes			
15ml Falcon tubes			
Incubator (37°C & 55°C)			
Beckman Centrifuge			
Amicon Ultra centrifuge tubes (15mL)	Fisher	UFC801096	
Phenol:Chloroform:IAA (25:24:1)	Sigma	77617-500ML	
Chloroform:IAA (24:1)	Sigma	25666-500ML	
TE buffer, pH 8.0			
1.5µl eppendorf tubes			

**Figures:**



**Figure 1.** 0.8% agarose gel electrophoresis image of high molecular weight DNA collected from four depths in the subarctic Pacific Ocean (in duplicate), stained with the intercalating agent ethidium bromine (10 mg/ml). Gel was run at 15V for ~16hrs in 1X TAE gel running buffer. Sample bands are of good quality showing little evidence of mechanical shearing (shows single bands or smears as opposed to multiple bands) although the 10m extracts retain some RNA carry over (see smear in the 0.5 to 2.0 Kb range).

## 16S rDNA Pyrotag Sequencing Protocol using KAPA2G Robust HotStart DNA Polymerase (Carmen Li, University of Alberta)

### Materials:

- 2.5  $\mu$ M stocks of primers (made up with sterile-filtered nuclease-free water):
- DMSO
- KAPA2G Robust HotStart DNA polymerase (5 U/ $\mu$ l)
- 5X KAPA2G Buffer A (included with KAPA2G *Taq* Polymerase)
- 5X KAPA Enhancer 1 (included with KAPA2G *Taq* Polymerase)
- 25mM MgCl<sub>2</sub> (included with KAPA2G *Taq* Polymerase)
- dNTP mix, 10 mM each (included with KAPA2G *Taq* Polymerase)
- Sterile nuclease-free water (Ambion)
- A thermocycler

### Method:

- 1) For each unique barcode (i.e., each different sample), add the following components together in an eppendorf tube to form a master mix:

	Per reaction:	Per 4 reactions:
2.5 $\mu$ M Forward primer	2.5 $\mu$ l	10 $\mu$ l
2.5 $\mu$ M Reverse primer	2.5 $\mu$ l	10 $\mu$ l
5X KAPA2G buffer A	5 $\mu$ l	20 $\mu$ l
(with MgCl <sub>2</sub> )		
5X KAPA Enhancer 1	5 $\mu$ l	20 $\mu$ l
DMSO	1.25 $\mu$ l	5 $\mu$ l
DNA polymerase	0.1 $\mu$ l	0.4 $\mu$ l
dNTP mix (10mM each)	0.5 $\mu$ l	2 $\mu$ l
MgCl <sub>2</sub> (25mM)	2 $\mu$ l	8 $\mu$ l
<u>Sterile nuclease-free water</u>	<u>5.15 <math>\mu</math>l</u>	<u>46 <math>\mu</math>l</u>
Total	24.5 $\mu$ l	98 $\mu$ l

- 2) Vortex each master mix. Dispense into 4 PCR tubes, one of these will be the no-template negative control.
- 3) Add 0.5  $\mu$ l of the same template to each of the 3 other PCR tubes.
- 4) Place all tubes in the thermocycler and run the following program:

1 cycle of 5 minutes denaturation at 95°C

10 cycles of 30seconds at 95°C, 30 seconds at 60°C, decreasing 0.5°C/cycle, and 30 seconds at 72°C

30 cycles of 30 seconds at 95°C, 30 seconds at 55°C and 30 seconds at 72°C.

1 cycle of 5 minutes at 72°C

- 5) Pool replicate template-containing PCR tubes. Run 4  $\mu$ l each template-containing and no-template control reaction on an agarose gel.
- 6) Provided bands are visible following agarose gel electrophoresis, purify the template-containing PCR reaction using the Qiagen QIAquick spin kit.
  - *I use the following modifications to increase template yield: Heat buffer EB to  $\sim 65^{\circ}\text{C}$  prior to elution; elute PCR products in 30  $\mu$ l buffer EB; following application of buffer EB to the column let it sit for 2 min before centrifugation.*
- 7) Quantify purified PCR products using Nanodrop. Perform at least 2 readings for each sample, mixing samples well before measuring.
  - *I take more than 2 readings if the first 2 differ by  $>2\text{ ng}/\mu\text{l}$ . If all replicates are relatively close I average all readings. If replicates are radically different (e.g., possibly from a different sample altogether) I omit the erroneous replicate.*
- 8) In 8-chain tubes or a 96-well plate, dilute the purified PCR products with PCR-grade water to the concentration of 30  $\text{ng}/\mu\text{l}$  or that of the most dilute sample in the sample set. I use 25  $\mu$ l final volumes.
- 9) Run 2  $\mu$ l of each diluted PCR product on an agarose gel.
- 10) Quantify diluted PCR products using Nanodrop.
- 11) After filling out the appropriate documentation, send the diluted PCR products for pyrotag sequencing.

## Genomic DNA Extraction from Tailing Ponds

### Overview

This protocol is a slight variation from methodology developed by Sangwon et al. (2009) and tailored to suit the need of obtaining high molecular weight DNA from low biomass environmental samples such as tailing ponds. The resulting DNA from this modified protocol has been shown to have a distribution of high molecular weight DNA essential for a construction of fosmid metagenomic library.

**Part I** describes the cell lysis step and the isolation of genomic DNA from the environmental matrix and cell debris. The cell lysis is achieved mainly through chemical lysis. The genomic DNA is then isolated through chloroform isoamyl alcohol and buffer exchanged using low molecular cut off filter. Mechanical lysis and aggressive vortexing are avoided in this procedure. This in combination with guanidine isothiocyanate and hexadecyltrimethylammonium bromide (CTAB), helps preserve the integrity of the high molecular weight genomic DNA. **Part II** describes the quality control step, where the resulting DNA is assessed for its quantity, quality, and purity. This would determine whether the DNA is appropriate for fosmid library construction.

### Part I: DNA Extraction

Time required: 6 hours

Before you start the extraction, pre-chill 50 ml tubes containing 20 ml chloroform-isoamyl alcohol (24:1) on ice. Set the oven to 65 °C to pre-heat. Check CTAB solution, if it is crystallized warm to 65°C to melt the crystals. Complete the denaturing solution and extraction buffer by adding the last ingredients just before you start the extraction (see the recipe).

1. Prepare the denaturing solution and extraction buffer.

Keep the denaturing solution on ice and the extraction buffer at 65°C (hybridization oven).

2. Aliquot 6 g of the tailing ponds material into 50 ml Falcon tube.
3. Add 3 ml of the denaturing solution (1 ml per 2 g of material).
4. Shake the extraction buffer before use. Add 15 ml of the extraction buffer. Shake briefly to mix. Store the remainder of the extraction buffer at 65°C.

Shaking the extraction buffer ensures homogeneity of the solution. Shaking the sample provides a minimum mechanical lysis and ensures most of the material is well permeated throughout the sample.

5. Incubate for 40 minutes at 65°C in a hybridization oven. During the incubation gently invert the tubes every 10 min, or continuously rotate the tubes at the lowest speed.

The heat helps weaken the cell membranes and aids in the cell lysis. The slight agitation help ensures the tailing material is well in contact with the extraction buffer.

6. Centrifuge at 1800 x g for 10 minutes using a swinging bucket rotor at 10 - 14°C.

JS 5.3 rotor from Beckman Coulter was used in this application.

7. Transfer the supernatant into a pre-chilled 50 ml Falcon tube containing 20 ml chloroform:IAA (24:1). Avoid transferring the detergent precipitate and the tailings material. Keep on ice.

When transferring the supernatant, it is better to leave some behind to ensure clean DNA.

8. Repeat the extraction process with another 15 ml of extraction buffer. Incubate at 65°C as previously described but only for 10 minutes. Centrifuge at 1800 x g for 10 minutes at 10 - 14°C.

Prior to the incubation, mix the pelleted tailings material with an autoclaved spatula for a better mixing.

9. Transfer the supernatant into the same tube as in step 7.

10. Cap the tube firmly and invert the tube to mix several times.

This is performed to increase the efficiency of the extraction by increasing the surface contact between the aqueous layer and the organic layer.

11. Place the tube over a bed of ice and gently shake the tube containing the supernatant chloroform:IAA mix using a rocking platform at low setting (1 - 2 out of 10) for 10 minutes. The agitation would provide constant mixing of the two phases further increasing the efficiency of the extraction. The low temperature would precipitate the detergent used and make it easier to be extracted by the organic phase.

12. Centrifuge at 1800 x g for 20 minutes at 10 - 14°C.

The spin would partition the phase and concentrating the cell debris in the interphase.

13. Transfer the aqueous phase (top) to a new 50 ml Falcon tube.

The organic extraction can be repeated to ensure clean DNA. However, multiple round of organic extraction would lower the final DNA yield.

When transferring the supernatant, it is better to leave some behind to ensure clean DNA. Make sure that the newly transferred aqueous phase is chloroform free. Residual chloroform would interfere with downstream steps.

14. Concentrate the aqueous phase using 10 kDa centrifugal filter unit (ie. Amicon filter).

In this application Amicon-15 with 10 kDa cut off was used. The spin time would vary across different sample from 10 to 20 minutes with a full load. The flow-through can be saved until the concentrate has been confirmed to contain DNA. The spin was done using the JS 5.3 swinging bucket rotor at 3000 xg for 15 minutes.

15. Buffer exchange the concentrate using 1x Tris EDTA buffer via diafiltration (top up the filter by topping up the filter and diluting the concentrate 1000x).

Buffer exchange is done to dilute the high salinity extraction buffer with storage buffer.

16. Transfer the concentrate into a microcentrifuge tube and note the volume.

## **Part II: Quality Control**

1. Quantify the DNA concentration via Quant-iT™ fluorometer.
2. Check for amplifiability through PCR reaction.
3. Run pulse field gel electrophoresis to check for the fragment distribution.

For Fosmid DNA library construction, the median size of the DNA should be equal to or exceed 36 kb.

## Supplementary Material

### \*Denaturing solution

4 M Guanidine isothiocyanate  
10 mM Tris-HCl (pH 7.0)  
1 mM EDTA  
0.5% 2-mercaptoethanol

### \*Note:

The denaturing solution is water based. Autoclave the solution excluding 2-mercaptoethanol and keep it at 4°C. The 2-mercaptoethanol is added just before use. If possible use fresh denaturing solution for the extraction, and the solution should not be stored for more than a week.

### \*\*Extraction Buffer

100 mM Sodium phosphate buffer [pH 7.0]  
100 mM Tris-HCl [pH 7.0]  
100 mM EDTA [pH 8.0]  
1.5 M NaCl  
1% Hexadecyltrimethylammonium bromide (CTAB)  
2% SDS

### \*\*Note:

The extraction buffer is water based. Autoclave and store the three solutions (salt buffer solution, CTAB, and SDS stock) separately. These solutions can be kept at room temperature. Combine the three components just before use. The precipitated detergent can be re-solubilized by heating them at 65°C.

Equipment or Reagent	Company	Catalogue Number	Comments
Centrifuge	Beckman Coulter	369003	Any centrifuge with 16000 ´ g capacity
Centrifuge rotor, JS5.3	Beckman Coulter	368690	For 1800 – 3000 g
Hybridization oven	Fisher Scientific	13-247-10	65 °C. It can be replaced with water bath.
Rocking plate	VWR	Rocking Platform Model 200	
1.5 ml tube	Eppendorf	0030 125.150	1 per sample
50 ml disposable tube	Corning	430290	3 per sample
Graduated glass cylinder			It can be replaced with glass pipette
Pipettes and filtered tips			20 ml, 200 ml, 1ml, 5 ml
Spatula			1 per sample, autoclaved
Chloroform-isoamyl alcohol (24:1)	Fluka	UN1888	20 ml per sample

Denaturing solution*			3 ml per sample
Guanidine isothiocyanate	Promega	V2791	
Trizma® base	Sigma	T1503	
Ethylenediaminetetraacetic acid disodium salt dihydrate	Sigma	E5134	
Mercaptoethanol	Sigma	M3148	
Extraction buffer**			30 ml per sample
Sodium phosphate dibasic	Sigma	S3264	
Sodium phosphate monobasic	Sigma	S3139	
Sodium Chloride	Sigma	S3014	
Hexadecyltrimethylammonium bromide	Sigma	H6269	
Quant-iT™ PicoGreen® dsDNA Reagent and Kit	Invitrogen	P758	
Spectrophotometer	Thermo Scientific	Nanodrop	
CHEF Mapper® XA Pulse Field Electrophoresis System	Bio-Rad	170-3670 to 170-3673	

### Reference:

Hurt, A. H., *et al.* Simultaneous Recovery of RNA and DNA from soils and sediments. AEM 67(10): 4495-4503 (2001).

Lee, S., Hallam, S. J., Extraction of High Molecular Weight Genomic DNA from Soils and Sediments. <http://www.jove.com/details.php?id=1569> doi: 10.3791/1569. J Vis Exp. 33 (2009))